

# FAS/CD95 antibody



Catalog Number: 107204

## Product name

FAS/CD95 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FAS/CD95 Mouse Monoclonal antibody. Positive FC detected in Jurkat cells. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of FAS/CD95 recombinant protein (Accession Number: NM\_000043). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107204(FAS antibody) at dilution of 1:500



$1 \times 10^6$  Jurkat cells were stained with 0.2ug FAS antibody (Catalog No:107204, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.