

# Fetuin-A antibody



Catalog Number: 107206

## Product name

Fetuin-A antibody

WB: 1:2000-1:20000

IP: 1:1000-1:10000

## Specificity

Human; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

Fetuin-A Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human liver cancer tissue, human liver tissue. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, human blood tissue, human heart tissue. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of Fetuin-A recombinant protein (Accession Number: NM\_001622). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

ELISA, WB, IF, IP, IHC


## Dilutions

Recommended Dilution:

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107206(AHSG antibody) at dilution of 1:5000




IP Result of anti-AHSG (IP:Catalog No:107206, 3ug; Detection:Catalog No:107206 1:3000) with HepG2 cells lysate 6000ug.


# Fetuin-A antibody

Catalog Number: 107206


---



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107206(AHSG antibody) at dilution of 1:200 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using AHSG antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107206(AHSG antibody) at dilution of 1:200 (under 10x lens)