

Fetuin-A antibody

Catalog Number: 107206

Product name

Fetuin-A antibody

WB: 1:2000-1:20000

IP: 1:1000-1:10000

Specificity

Human; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

Fetuin-A Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human liver cancer tissue, human liver tissue. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, human blood tissue, human heart tissue. Observed molecular weight by Western-blot: 58 kDa

Validations



Preparation

This antibody was obtained by immunization of Fetuin-A recombinant protein (Accession Number: NM_001622). Purification method: Protein A purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107206(AHSG antibody) at dilution of 1:5000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type


Mouse IgG2a

Applications

ELISA, WB, IF, IP, IHC

Dilutions


Recommended Dilution:




IP Result of anti-AHSG (IP:Catalog No:107206, 3ug; Detection:Catalog No:107206 1:3000) with HepG2 cells lysate 6000ug.

Fetuin-A antibody

Catalog Number: 107206



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107206(AHSG antibody) at dilution of 1:200 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using AHSG antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107206(AHSG antibody) at dilution of 1:200 (under 10x lens)

