

EEF1D antibody



Catalog Number: 107249

Product name

EEF1D antibody

Specificity

human, mouse; other species not tested.

Antibody description

EEF1D Mouse Monoclonal antibody. Positive IF detected in MCF-7. Positive IHC detected in human pancreas cancer. Observed molecular weight by Western-blot: 31/71kd

Preparation

This antibody was obtained by immunization of EEF1D recombinant protein (Accession Number: XM_047421433). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:107249(EEF1D antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:107249(EEF1D antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells, using EEF1D antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).