

HEXA antibody



Catalog Number: 107261

Product name

HEXA antibody

Specificity

Human; other species not tested.

Antibody description

HEXA Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of HEXA recombinant protein (Accession Number: NM_000520). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107261(HEXA Antibody) at dilution of 1:16000