

Gastrin antibody

Catalog Number: 107291

Product name

Gastrin antibody

Specificity

Human; other species not tested.

Antibody description

Gastrin Mouse Monoclonal antibody. Positive IHC detected in human stomach tissue.

Preparation

This antibody was obtained by immunization of Gastrin recombinant protein (Accession Number: NM_000805). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:100-1:700

Validations



Immunohistochemistry of paraffin-embedded human stomach tissue slide using (GAST Antibody) at dilution of 1:1000 (under 10x lens).



Immunohistochemistry of paraffin-embedded human stomach tissue slide using (GAST Antibody) at dilution of 1:1000 (under 40x lens).