

# NCALD antibody



Catalog Number: 107323

## Product name

NCALD antibody

## Specificity

Human; other species not tested.

## Antibody description

NCALD Mouse Monoclonal antibody. Positive WB detected in Human brain, human brain tissue. Positive IHC detected in human cerebellum tissue, human testis tissue. Observed molecular weight by Western-blot: 22kD

## Preparation

This antibody was obtained by immunization of NCALD recombinant protein (Accession Number: XM\_017013902). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107323(NCALD antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human cerebellum using Catalog No:107323(NCALD antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:107323(NCALD Antibody) at dilution of 1:50