

NDUFS5 antibody



Catalog Number: 107326

Product name

NDUFS5 antibody

Specificity

Human; other species not tested.

Antibody description

NDUFS5 Mouse Monoclonal antibody. Positive WB detected in Human kidney. Positive IHC detected in human kidney tissue, human liver tissue. Observed molecular weight by Western-blot: 15kD

Preparation

This antibody was obtained by immunization of NDUFS5 recombinant protein (Accession Number: NM_001184979). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107326(NDUFS5 antibody) at dilution of



1:1000 Immunohistochemical of paraffin-embedded human kidney using Catalog No:107326(NDUFS5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107326(NDUFS5 antibody) at dilution of 1:50 (under 40x lens)