

NOB1 antibody



Catalog Number: 107334

Product name

NOB1 antibody

Immunogen

[Human NOB1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

NOB1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, human placenta tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of NOB1 recombinant protein (Accession Number: BC000050). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations




Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107334(NOB1 Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107334(NOB1 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).

NOB1 antibody

Catalog Number: 107334



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107334(NOB1 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).