

NR1H2 antibody



Catalog Number: 107335

Product name

NR1H2 antibody

Specificity

Human, Pig; other species not tested.

Antibody description

NR1H2 Mouse Monoclonal antibody. Positive WB detected in human liver tissue, pig liver tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of NR1H2 recombinant protein (Accession Number: NM_007121). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107335(NR1H2 Antibody) at dilution of 1:1000