

# NR1H2 antibody



Catalog Number: 107335

---

## Product name

NR1H2 antibody

## Specificity

Human, Pig; other species not tested.

## Antibody description

NR1H2 Mouse Monoclonal antibody. Positive WB detected in human liver tissue, pig liver tissue. Observed molecular weight by Western-blot: 51 kDa

## Preparation

This antibody was obtained by immunization of NR1H2 recombinant protein (Accession Number: NM\_007121). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107335(NR1H2 Antibody) at dilution of 1:1000