

NR1H3 antibody



Catalog Number: 107336

Product name

NR1H3 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

NR1H3 Mouse Monoclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in MCF7 cells, HepG2 cells. Observed molecular weight by Western-blot: 50 kDa 44 kDa

Preparation

This antibody was obtained by immunization of NR1H3 recombinant protein (Accession Number: XM_047426198). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107336(NR1H3 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:107336(NR1H3 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)