

NR2F6 antibody



Catalog Number: 107338

Product name

NR2F6 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NR2F6 Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of NR2F6 recombinant protein (Accession Number: NM_005234). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107338(NR2F6 antibody) at dilution of 1:50 (under 40x lens)



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107338(NR2F6 Antibody) at dilution of 1:2000