

IMP3 antibody



Catalog Number: 107347

Product name

IMP3 antibody

Specificity

Human; other species not tested.

Antibody description

IMP3 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, SKOV-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of IMP3 recombinant protein (Accession Number: BC020557). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107347(IMP3 Antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells using Catalog No:107347(IMP3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)