

Inhibin beta A antibody



Catalog Number: 107349

Product name

Inhibin beta A antibody

Specificity

Human; other species not tested.

Antibody description

Inhibin beta A Mouse Monoclonal antibody.
Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Inhibin beta A recombinant protein (Accession Number: XM_017012176). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107349(INHBA antibody) at dilution of 1:500