

# Inhibin beta A antibody



Catalog Number: 107349

## Product name

Inhibin beta A antibody

## Specificity

Human; other species not tested.

## Antibody description

Inhibin beta A Mouse Monoclonal antibody.  
Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of Inhibin beta A recombinant protein (Accession Number: XM\_017012176). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107349(INHBA antibody) at dilution of 1:500