

# KCNN4 antibody



Catalog Number: 107360

## Product name

KCNN4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KCNN4 Mouse Monoclonal antibody. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of KCNN4 recombinant protein (Accession Number: NM\_002250). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107360(KCNN4 Antibody) at dilution of 1:1000