

KIFAP3 antibody



Catalog Number: 107362

Product name

KIFAP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KIFAP3 Mouse Monoclonal antibody. Positive WB detected in human testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 91-100kd

Preparation

This antibody was obtained by immunization of KIFAP3 recombinant protein (Accession Number: NM_014970). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107362(KIFAP3 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:107362(KIFAP3 Antibody) at dilution of 1:50