

L2HGDH antibody



Catalog Number: 107368

Product name

L2HGDH antibody

Immunogen

[Human L2HGDH Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

L2HGDH Mouse Monoclonal antibody. Positive WB detected in pig brain tissue, MCF-7 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of L2HGDH recombinant protein (Accession Number: BC006117). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



pig brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107368(L2HGDH Antibody) at dilution of 1:1000