

LAMP2 antibody



Catalog Number: 107370

Product name

LAMP2 antibody

Specificity

Human; other species not tested.

Antibody description

LAMP2 Mouse Monoclonal antibody. Positive WB detected in HUVEC cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of LAMP2 recombinant protein (Accession Number: NM_013995). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgA

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:107370(LAMP2 Antibody) at dilution of 1:1000