

IL12A antibody



Catalog Number: 107379

Product name

IL12A antibody

Specificity

Human; other species not tested.

Antibody description

IL12A Mouse Monoclonal antibody. Positive WB detected in Recombinant protein.

Preparation

This antibody was obtained by immunization of IL12A recombinant protein (Accession Number: NM_000882). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107379(IL12A Antibody) at dilution of 1:2000