

IL-17D antibody

Catalog Number: 107384

Product name

IL-17D antibody

Specificity

Human; other species not tested.

Antibody description

IL-17D Mouse Monoclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of IL-17D recombinant protein (Accession Number: NM_138284). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107384(IL17D antibody) at dilution of 1:20000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107384(IL17D Antibody) at dilution of 1:1000