

# IL-17D antibody



Catalog Number: 107384

## Product name

IL-17D antibody

## Specificity

Human; other species not tested.

## Antibody description

IL-17D Mouse Monoclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of IL-17D recombinant protein (Accession Number: NM\_138284). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107384(IL17D antibody) at dilution of 1:20000



	-	+	-	+	
	-	+	-	+	Overexpression of IL17D his-myc
	+	+	-	-	Mouse Anti human IL17D monoclonal antibody
	-	-	+	+	Mouse Anti his-tag monoclonal antibody

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107384(IL17D Antibody) at dilution of 1:1000