

# IL1R2 antibody

Catalog Number: 107389

## Product name

IL1R2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL1R2 Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 45KD

## Preparation

This antibody was obtained by immunization of IL1R2 recombinant protein (Accession Number: NM\_173343). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107389(IL1R2 antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:107389(IL1R2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)