

# IL21R antibody



Catalog Number: 107392

## Product name

IL21R antibody

## Specificity

Human; other species not tested.

## Antibody description

IL21R Mouse Monoclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in K-562 cells, Jurkat cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of IL21R recombinant protein (Accession Number: NM\_181078). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:200


## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107392(IL21R Antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (IL21R Antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using (IL21R Antibody) at dilution of 1:100 (under 40x lens)