

# IL22RA2 antibody

Catalog Number: 107393

## Product name

IL22RA2 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL22RA2 Mouse Monoclonal antibody. Positive IHC detected in human small intestine tissue, human breast cancer tissue, human lung cancer tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of IL22RA2 recombinant protein (Accession Number: NM\_181309). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107393(IL22RA2 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:107393(IL22RA2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog

# IL22RA2 antibody



Catalog Number: 107393

---

No:107393(IL22RA2 Antibody) at dilution of 1:50 (under 40x lens)