

IL-27 antibody

Catalog Number: 107395

Product name

IL-27 antibody

Specificity

Human; other species not tested.

Antibody description

IL-27 Mouse Monoclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 31 kDa

Preparation

This antibody was obtained by immunization of IL-27 recombinant protein (Accession Number: NM_145659). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

IHC: 1:20-1:200


Validations



IL27 fusion protein were subjected to SDS PAGE followed by western blot with Catalog No:107395(IL27 Antibody) at dilution of 1:5000




Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107395(IL27 Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107395(IL27 Antibody) at dilution of 1:50 (under 10x lens)

IL-27 antibody

Catalog Number: 107395



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107395(IL27 Antibody) at dilution of 1:50 (under 40x lens)