

# IL-28A antibody

Catalog Number: 107396

## Product name

IL-28A antibody

## Specificity

Human; other species not tested.

## Antibody description

IL-28A Mouse Monoclonal antibody. Positive WB detected in ag18525. Positive IHC detected in human tonsillitis tissue, human small intestine tissue. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of IL-28A recombinant protein (Accession Number: NM\_172138). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:5000-1:50000

IHC: 1:20-1:200


## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107396(IL28A Antibody) at dilution of 1:5000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107396(IL28A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107396(IL28A Antibody) at dilution of 1:50

# IL-28A antibody



Catalog Number: 107396

---

(under 40x lens)