

# IL-29 antibody



Catalog Number: 107397

## Product name

IL-29 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL-29 Mouse Monoclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of IL-29 recombinant protein (Accession Number: NM\_172140). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:5000-1:50000

## Validations



GST tagged recombinant IL29 protein were subjected to SDS PAGE followed by western blot with Catalog No:107397(IL29 Antibody) at dilution of 1:20000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107397(IL29 Antibody) at dilution of 1:5000