

LIPH antibody



Catalog Number: 107405

Product name

LIPH antibody

Specificity

Human; other species not tested.

Antibody description

LIPH Mouse Monoclonal antibody. Positive WB detected in SW 1990 cells, BxPC-3 cells, COLO 320 cells. Positive IHC detected in human pancreas tissue, human colon cancer tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of LIPH recombinant protein (Accession Number: NM_139248). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:50-1:500

Validations



SW 1990 cells were subjected to SDS PAGE followed by western blot with Catalog No:107405(LIPH Antibody) at dilution of 1:2000