

LPCAT1 antibody

Catalog Number: 107408

Product name

LPCAT1 antibody

Specificity

human; other species not tested.

Antibody description

LPCAT1 Mouse Monoclonal antibody. Positive WB detected in MCF7 cells. Positive IF detected in HeLa. Observed molecular weight by Western-blot: 55kD

Preparation

This antibody was obtained by immunization of LPCAT1 recombinant protein (Accession Number: BC020166). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107408(LPCAT1 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells, using LPCAT1 antibody Catalog No: at 1:50 dilution and Rhodamine-labeled goat anti-mouse IgG (red).