

# MLANA antibody

Catalog Number: 107426

## Product name

MLANA antibody

## Specificity

Human; other species not tested.

## Antibody description

MLANA Mouse Monoclonal antibody. Positive IHC detected in human malignant melanoma tissue.

## Preparation

This antibody was obtained by immunization of MLANA recombinant protein (Accession Number: NM\_005511). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human malignant melanoma tissue slide using Catalog No:107426(MLANA Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human malignant melanoma tissue slide using Catalog No:107426(MLANA Antibody) at dilution of 1:200 (under 40x lens).