

# PA2G4 antibody



Catalog Number: 107448

## Product name

PA2G4 antibody

## Specificity

Human; other species not tested.

## Antibody description

PA2G4 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, MCF7 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 38-42kD, 48kD

## Preparation

This antibody was obtained by immunization of PA2G4 recombinant protein (Accession Number: NM\_006191). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107448(PA2G4 antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells, using PA2G4 antibody Catalog No: at 1:50 dilution and Rhodamine-labeled goat anti-mouse IgG (red).