

PAI-1 antibody

Catalog Number: 107449

Product name

PAI-1 antibody

Specificity

Human, Pig; other species not tested.

Antibody description

PAI-1 Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells, human plasma tissue, pig liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of PAI-1 recombinant protein (Accession Number: NM_000602). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107449(SERPINE1 Antibody) at dilution of 1:2000



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:107449(SERPINE1 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)