

PAX1-Specific antibody



Catalog Number: 107452

Product name

PAX1-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PAX1-Specific Mouse Monoclonal antibody.
Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_006192).
Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107452(PAX1-Specific Antibody) at dilution of 1:1000