

# PAX4 antibody

Catalog Number: 107453

## Product name

PAX4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PAX4 Mouse Monoclonal antibody. Positive WB detected in Hela cells, K-562 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 38 kDa, 58 kDa

## Preparation

This antibody was obtained by immunization of PAX4 recombinant protein (Accession Number: NM\_006193). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107453(PAX4 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:107453(PAX4 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)