

phospho(409/410)-TDP43 antibody



Catalog Number: 107470

Product name

Monoclonal

phospho(409/410)-TDP43 antibody

Ig Type

Mouse IgG1

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

phospho(409/410)-TDP43 Mouse Monoclonal antibody. Positive WB detected in Ethacrynic acid treated HeLa cells. Observed molecular weight by Western-blot: 45-50, 35 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007375).
Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Untreated and ethacrynic acid treated HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107470 (phospho(409/410)-TDP43 Antibody) at dilution of 1:850

Clonality