

# PHOX2B antibody



Catalog Number: 107471

## Product name

PHOX2B antibody

## Immunogen

[Human PHOX2B Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

PHOX2B Mouse Monoclonal antibody. Positive WB detected in Neuro-2a cells. Positive IHC detected in neuroblastoma tissue. Observed molecular weight by Western-blot: 32-36 kDa

## Preparation

This antibody was obtained by immunization of PHOX2B recombinant protein (Accession Number: NM\_003924). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:107471(PHOX2B Antibody) at dilution of



Immunohistochemistry of paraffin-embedded neuroblastoma tissue slide using (PHOX2B Antibody) at dilution of 1:500 (under 40x lens).