

# PPARD antibody



Catalog Number: 107480

---

## Product name

PPARD antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PPARD Mouse Monoclonal antibody. Positive WB detected in Colo320 cells, human heart tissue. Observed molecular weight by Western-blot: 54kd

## Preparation

This antibody was obtained by immunization of PPARD recombinant protein (Accession Number: NM\_177435). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:107480(PPARD antibody) at dilution of 1:1000