

# Protein C inhibitor antibody



Catalog Number: 107487

## Product name

Protein C inhibitor antibody

## Specificity

Human; other species not tested.

## Antibody description

Protein C inhibitor Mouse Monoclonal antibody. Positive WB detected in Human testis, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of Protein C inhibitor recombinant protein (Accession Number: NM\_000624). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107487(SERPINA5 antibody) at dilution of 1:1000