

Protein C inhibitor antibody



Catalog Number: 107487

Product name

Protein C inhibitor antibody

Specificity

Human; other species not tested.

Antibody description

Protein C inhibitor Mouse Monoclonal antibody. Positive WB detected in Human testis, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of Protein C inhibitor recombinant protein (Accession Number: NM_000624). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107487(SERPINA5 antibody) at dilution of 1:1000