

# Rab18 antibody



Catalog Number: 107491

## Product name

Rab18 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Rab18 Mouse Monoclonal antibody. Positive WB detected in Human adrenal gland, mouse testis tissue. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of Rab18 recombinant protein (Accession Number: NM\_021252). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:107491(RAB18 antibody) at dilution of 1:1000