

RABEP2 antibody

Catalog Number: 107493

Product name

RABEP2 antibody

Specificity

Human; other species not tested.

Antibody description

RABEP2 Mouse Monoclonal antibody. Positive WB detected in PC-3 cells. Positive IHC detected in human prostate cancer tissue, human breast cancer tissue. Observed molecular weight by Western-blot: 62kd

Preparation

This antibody was obtained by immunization of RABEP2 recombinant protein (Accession Number: BC058900). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107493(RABEP2 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:107493(RABEP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog

RABEP2 antibody



Catalog Number: 107493

No:107493(RABEP2 Antibody) at dilution of 1:50 (under 40x lens)