

# Rac1 antibody



Catalog Number: 107494

## Product name

Rac1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Rac1 Mouse Monoclonal antibody. Positive IP detected in HL-60 cells. Positive WB detected in HL-60 cells, human brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of Rac1 recombinant protein (Accession Number: NM\_006908). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:107494(RAC1 antibody) at dilution of 1:1000



IP Result of anti-RAC1 (IP:Catalog No:107494, 3ug; Detection:Catalog No:107494 1:700) with HL-60 cells lysate 4000ug.