

RAC2 antibody

Catalog Number: 107495

Product name

RAC2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RAC2 Mouse Monoclonal antibody. Positive IHC detected in human spleen tissue, human tonsillitis tissue. Positive WB detected in HL60 cells, Jurkat cells. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of RAC2 recombinant protein (Accession Number: NM_002872). Purification method: Protein G purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:107495(RAC2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human spleen slide using Catalog No:107495(RAC2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human spleen slide using Catalog No:107495(RAC2 Antibody) at dilution of 1:50