

RBBP9 antibody

Catalog Number: 107499

Product name

RBBP9 antibody

Specificity

Human; other species not tested.

Antibody description

RBBP9 Mouse Monoclonal antibody. Positive WB detected in A549 cell. Positive IHC detected in human pancreas tissue, human pancreas cancer tissue. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of RBBP9 recombinant protein (Accession Number: NM_006606). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:107499(RBBP9 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:107499(RBBP9 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:107499(RBBP9 antibody) at dilution of 1:50 (under 40x lens)