

# RNH1 antibody



Catalog Number: 107508

---

## Product name

RNH1 antibody

## Immunogen

[Human RNH1 Recombinant protein \(His tag\)](#)

## Specificity

human; other species not tested.

## Antibody description

RNH1 Mouse Monoclonal antibody. Positive WB detected in lo2 cells. Observed molecular weight by Western-blot: 50kD

## Preparation

This antibody was obtained by immunization of RNH1 recombinant protein (Accession Number: NM\_203388). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107508(RNH1 antibody) at dilution of 1:500