

RXRA antibody



Catalog Number: 107514

Product name

Monoclonal

RXRA antibody

Ig Type

Specificity

Mouse IgG1

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

RXRA Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 44 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of RXRA recombinant protein (Accession Number: BC007925). Purification method: Protein G purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107514(RXRA antibody) at dilution of 1:100