

# RXRA antibody



Catalog Number: 107514

---

## Product name

RXRA antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RXRA Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 44 kDa

## Preparation

This antibody was obtained by immunization of RXRA recombinant protein (Accession Number: BC007925). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107514(RXRA antibody) at dilution of 1:100