

SALL4 antibody



Catalog Number: 107520

Product name

SALL4 antibody

Specificity

Human; other species not tested.

Antibody description

SALL4 Mouse Monoclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 165 kDa

Preparation

This antibody was obtained by immunization of SALL4 recombinant protein (Accession Number: NM_020436). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107520(SALL4 Antibody) at dilution of 1:1000